

TAPXORA LIMS APP

Detailed Training Guide and Demo Deck

The Whole Journey In One View


- 1. User signs in
- 2. Patient is registered
- 3. Order is created with one or more tests
- 4. Sample label is generated and attached
- 5. Sample is received, stored, and tracked
- 6. Result is entered by laboratory staff
- 7. Result is verified by an authorized verifier
- 8. Report is printed or downloaded
- 9. Invoice and receipt are handled
- 10. Dashboard, stock, and audit logs stay updated

Training Tip


The easiest demo is to follow one patient from account access, to patient registration, to order creation, to sample tracking, to result entry, to report printing, and finally to billing.

Entry Screens For New Users

These are real screenshots from the live app. They are useful during training because they show the exact first pages staff will see when they open the system.



Login page: users can sign in with email and password or request a magic link.



Register page: new staff enter their name, email, and password to create an account.

User Roles And What Each Role Does

A good training session should explain that users do not all see the same tools. The system shows menus based on responsibility.

Admin

Creates tests, assigns roles, reviews audit logs, manages settings, and oversees the whole laboratory workflow.

Receptionist

Creates patient records, opens lab orders, confirms patient details, and supports front desk registration.

Lab Scientist

Receives samples, updates sample status, enters test results, and supports stock usage logging.

Verifier

Checks entered results, confirms they are safe to release, and approves final reports.

Accountant

Monitors invoices, records payments, prints receipts, and reviews daily revenue summaries.

1. Introduction: What TAPXORA LIMS APP Does

This app helps a laboratory run its work from beginning to end. It keeps patient records, tracks samples, saves results, prints reports, records payments, and watches inventory in one connected system.

Step by step

- Think of the app as the lab's organized notebook, calculator, tracker, and report printer combined into one safe digital workspace.
- Every action in the lab moves through a clear path: patient, order, sample, result, report, and payment.
- Because the app is offline first, staff can keep working when the internet is poor and sync later.

Helpful notes

- Use this guide during staff onboarding, product demonstrations, and refresher training.
- The easiest way to teach the app is to follow one sample patient through the full workflow.

2. User Registration and Sign In

Nobody can use the app until they have an account. The entry screens are the first place new users learn how to enter the system.

Step by step

- Open the Register page.
- Type the user's full name, email address, and password.
- Click Create Account to save the new login.
- The user can also request a magic link if passwordless sign in is preferred.
- After registration, the user can return to the Login page and sign in.

Helpful notes

- Use work email addresses so each account is easy to trace.
- If a role changes after the user is logged in, ask the user to sign out and sign in again.

3. Roles and Access Control

People in the laboratory do different jobs. The app uses roles so each person sees the tools that match their responsibilities.

Step by step

- Find the user in the profiles table or admin user management area.
- Assign one role: Admin, Receptionist, LabScientist, Verifier, or Accountant.
- Confirm the user belongs to the correct facility so records stay properly separated.
- Test the role by signing in and checking that the correct menu items appear.

Helpful notes

- Admins should be few because they can control major settings.
- Role based screens reduce confusion and lower the risk of accidental changes.

4. Patient Registration and Consent

A clean patient record makes every later step safer. If the patient identity is wrong, the sample and result workflow can also become wrong.

Step by step

- Open the Patients page and click New Patient.
- Enter the patient's name, phone number, sex, date of birth, and address.
- Add extra details such as emergency contact, notes, and local identifiers when available.
- Tick the NDPR consent checkbox before saving.
- Save the patient and confirm the patient appears in the search list.

Helpful notes

- Use the phone number carefully because it helps staff find returning patients quickly.
- Do not create duplicate patient records unless you have confirmed the person is truly new.

5. Patient Search and History

When a patient returns, staff should not start from zero. Search and history help the team find old records, repeat orders, and earlier reports.

Step by step

- Open the patient list.
- Search by patient name, phone number, or lab ID.
- Open the patient profile.
- Review previous orders, previous results, and earlier billing records if needed.

Helpful notes

- History is especially useful for follow up care and repeat investigations.
- Use filters when the patient list grows large.

6. Test Catalogue Management

The test catalogue is the master list of tests the lab offers. It controls pricing, result formats, and reference ranges.

Step by step

- Open the Test Catalogue module as an Admin.
- Create a test and enter the test name and price.
- Choose the result type such as numeric, text, dropdown, or positive or negative.
- Enter a reference range using min and max values or plain text when needed.
- Set the test to active so it becomes available during order creation.

Helpful notes

- Inactive tests stay hidden from new orders but remain in history.
- Correct pricing here helps billing remain automatic and accurate.

7. Creating a Lab Order

An order connects the patient to one or more requested tests. This is the official start of the laboratory job.

Step by step

- Open the Orders page.
- Search for and select the patient.
- Choose one or more tests from the active catalogue.
- Set priority such as routine, urgent, or stat if your workflow uses it.
- Save the order so sample items can be generated.

Helpful notes

- Double check the patient before saving.
- One patient can have multiple tests inside the same order.

8. Sample Collection, Barcode Labels, and Storage

Samples must be labelled correctly and stored correctly. This reduces mistakes and supports chain of custody tracking.

Step by step

- After the order is created, print the generated sample barcode or QR label.
- Attach the label to the correct sample container immediately.
- Confirm the label matches the patient name, order, and requested test.
- Record where the sample is kept, such as bench, fridge, freezer, tray, rack, or shelf.
- Add collection notes when necessary, for example if the sample was difficult to collect or delayed.

Helpful notes

- Never place an unlabeled sample aside to label later.
- Good storage notes make it easier to locate the sample during follow up work.

9. Sample Reception, Tracking, and Chain of Custody

The lab should always know where the sample is and who handled it. This is important for trust, traceability, and quality control.

Step by step

- Open Sample Reception.
- Scan the barcode or type the sample code.
- Confirm the patient, order number, and test.
- Move the sample through status stages such as Registered, Collected, In Progress, Results Entered, Verified, and Reported.
- Allow the system to log each action in the custody history.

Helpful notes

- If a sample issue is noticed, record it immediately instead of waiting.
- Staff should update status as work happens so the dashboard stays useful.

10. Results Entry

This is where the scientific work is entered into the system. Clean result entry supports good reporting and good verification.

Step by step

- Open the Results Entry workspace.
- Choose the correct sample and test.
- Use the correct form type based on the test setup: numeric, text, dropdown, or positive or negative.
- Review any automatic abnormal flag generated from the reference range.
- Save the result for verification.

Helpful notes

- Always confirm the patient and sample before typing a result.
- Use comments when the result needs extra explanation.

11. Result Verification and Approval

Verification is a safety checkpoint. A second authorized person confirms the result before the report is released.

Step by step

- Open the verification queue.
- Review the entered value, comments, and abnormal flag.
- Compare the result to the correct patient and sample.
- Approve the result if everything is correct.
- If something is wrong, reject it and return it for correction.

Helpful notes

- This two step process improves quality and lowers reporting errors.
- Verified results become ready for professional report generation.

12. Report Printing and Download

Verified results can be turned into clean medical reports for clinicians and patients.

Step by step

- Open the Reports module.
- Choose a verified order.
- Review the patient details, test names, values, units, flags, and comments.
- Download the PDF report or use the print action.
- Release the report only after final review.

Helpful notes

- Printed reports should be easy to read and professionally formatted.
- Use the report view during demos because it helps users see the final value of the system.

13. Billing, Invoices, and Receipts

Billing keeps the business side of the lab connected to the service side. Orders, prices, invoices, and payments stay linked together.

Step by step

- Open the Billing module.
- Review the invoice automatically created from the ordered tests.
- Check the payment status: unpaid, partial, or paid.
- Record the amount received and payment method.
- Generate or print a receipt after payment is saved.

Helpful notes

- Automatic pricing works best when the test catalogue is accurate.
- Daily revenue summaries help supervisors review business performance.

14. Inventory and Reagent Tracking

The lab cannot work smoothly when reagents and consumables run out. Inventory makes shortages visible early.

Step by step

- Open Inventory and create items such as reagents, kits, tubes, gloves, and consumables.
- Enter quantity, lot number, expiry date, and reorder level.
- Record stock in when new items arrive.
- Record stock out or usage when staff use items during daily work.
- Watch low stock and near expiry alerts on the dashboard.

Helpful notes

- Lot numbers are useful for traceability during quality review.
- Do not wait until stock finishes before updating quantities.

15. Dashboards, Audit Logs, and Offline Work

Managers need a clear summary of operations, and the system needs to remain useful when internet quality is weak.

Step by step

- Use the dashboard to see today's worklist, turnaround time, test volume, revenue, and alerts.
- Use the audit logs viewer to see who changed what and when.
- When the internet drops, continue working with the offline mode.
- Allow queued changes to sync when connectivity returns.

Helpful notes

- Audit logs help during supervision, investigations, and compliance review.
- Offline support is especially important in real Nigerian lab environments.

16. Full Training Demo Script

The best way to teach the app is to act out a complete patient journey from start to finish.

Step by step

- Create one user account and explain what role that user has.
- Register one patient and explain the consent checkbox.
- Create one order with two tests.
- Show the barcode label and explain where the sample is stored.
- Move the sample through status updates.
- Enter one result, verify it, print the report, and record payment.
- End by showing inventory alerts, dashboard summaries, and the audit log.

Helpful notes

- This simple story based demo is easy for first time users to remember.
- Even a child can follow the workflow when each step is shown in order.

Simple End To End Demo Walkthrough

Use this page when you want to show the app quickly during a meeting, training class, or client presentation. Each card represents one simple action in the order people should learn it.

Step 1

Register a user account or sign in with a prepared account.

Step 2

Explain the user's role and which menu items that role can see.

Step 3

Register one patient and tick the NDPR consent checkbox.

Step 4

Create one lab order with one or two tests.

Step 5

Show the barcode label and explain sample storage location.

Step 6

Update sample status and explain chain of custody.

Step 7

Enter one result and point out abnormal flags.

Step 8

Verify the result, print the report, and record payment.